MODEL CA6 - HARDNESS ANALYZER

Compact online colorimeter for the automatic measurement of Hardness in water

APPLICATION FIELDS

- Power plants
- Cooling water
- Water steam cycle
- Boiler feedwater
- Reversed osmosis
- Ion exchangers
- Ultrapure water
- Drinking water



ADVANTAGES / FEATURES

Dual compartment enclosure

To ensure complete separation between the electronics and the wet part.

Low reagent consumption

Minimum operating cost by small reagent consumption, only 1.7L (0.45 US.gal) for the 16 mm cell / 2.5L (0.66 US.gal) for the 26 mm cell of each reagent every 30 days with 15 minute analysis frequency.

Automatic calibration / validation / cleaning

Validation, cleaning and calibration are standard features which significantly reduce downtime and operator intervention Multiple streams ensuring the most accurate results are obtained.

Free selectable validation, cleaning and calibration intervals.


Color touchscreen user interface

The CA6 Colorimeter is equipped with a graphic touchscreen interface showing measured values and status information. Easy access to menus and functions. Multiple languages. Integrated datalogger with USB download.

Factory tested, ready for installation and operation

Just connect the power, sample, and reagent lines and the analyzer is fully operational.

Dual streams version available. External Sequencer, switching up to 4 sample streams.



MEASUREMENT PRINCIPLE

The photometric determination is based on the reaction of calcium with O-cresolphalein Complexone solution, which yields a violet colored complex. The intensity of the color formed is proportional to the calcium concentration in the sample. Absorbance of the complex is measured at 572 nm.



TECHNICAL SPECIFICATIONS

Hardness as CaCO₃ (ppb, ppm, mg/l). Measured parameter:

Differential photometric absorbance. Measuring principle:

O-cresolphalein Complexone.

0-500 ppb (26 mm cell) - 0-1000 ppb (16 mm Measuring range:

cell), up to 50 ppm with internal dilution.

 \pm 5 ppb or \pm 5%, whichever is greater (26

mm cell)

Reproducibility: ± 10 ppb or ± 5%, whichever is greater (16

mm cell).

Freely programmable, batch near-continuous Analysis frequency:

analysis.

6 minutes, including conditioning before Cycle time:

analysis cycle and rinsing after measuring.

Temperature heated Reaction cell:

Pressure-free from overflow vessel

Temperature: 41 - 113 °F (5 to 45 °C) Sample:

Flow Rate: 80 to 500 mL/min Connection: 6 mm (1/4-in.)

Pressure-free, atmospheric drain Drain:

Connection: 12 mm (1/2-in.)

1, 2 with integrated switching valve N° of streams: 3, 4 with external sequencer

Dimensions (H x W x D): 23.6 x 15.0 x 8.2 in / 606 x 380 x 209 mm

Weight: Approx. 44 lbs (20 Kg)

Voltage: 100 - 240 VAC 50/60 Hz standard or 24 Power supply:

VDC (option)

Power consumption: max. 80 VA

2 x 4-20 mA outputs for measured data Outputs:

Modbus RTU RS485

Alarms: 4 SPDT programmable potential free relays

Digital input: Remote start / stop

41 - 113 °F (5 - 45 °C) Operating Temperature:

10 to 90% non-condensing (indoor use, Humidity:

outdoor installation only possible with protective cabinet or shelter not included)

Wall mount (standard), bench top support or

panel mount (options).

Ingress Protection: **IP54**

Installation:

